

Connect App Firmware Update Tool

- **Note:** Images shown are on iOS in this guide. Over-the-air (OTA) app is available on Google Play Store and App Store.



FIRST:



Ensure your mobile device has bluetooth turned on:
Go to your phone or tablet's **Settings > Bluetooth**, and turn Bluetooth **ON**.

- If the firmware app is not downloaded on your device, follow the process in Step 1. Otherwise proceed to Step 2.

Step 1:

Launch the Connect App for Quick Link to the Firmware Update Tool

Click the menu icon to open the navigation panel in the Connect App for a link to the applicable app store where installation is available.



Step 2:

Open OTA App

Long press the Bluetooth button on the Bluetooth capable device.




Step 3:

Connect


Select the serial number that appears on the app.

Tip: If no serial number appears and Bluetooth is displayed on the device screen, then the device does not need to be updated.



Step 4: Update the Device

If a firmware update is available, the information will display on the screen. Tapping update will start the process.



Update Successful!

